



THE GLOBE PROGRAM

A worldwide science and education program



Hydrosphere



Water pH Protocol
Using pH Paper





Overview Water pH Protocol using pH Paper

This module:

- Reviews the selection of a GLOBE hydrology site
- Reviews the water sampling technique used in GLOBE hydrology protocols
- Provides a step by step introduction of the protocol method

Learning Objectives

After completing this module, you will be able to:

- Define water pH and explain how environmental variables result in different measurements
- Describe the importance of instrument calibration in the the collection of accurate data
- Conduct water pH measurements using a pH meter
- Upload data to the GLOBE portal
- Visualize data using GLOBE's Visualization System

Estimated time needed to complete this module: 1.5 hours

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself


H. Additional resources



The Hydrosphere as part of the Earth System

The hydrosphere is the part of the Earth system that includes water, ice and water vapor. Water participates in many important natural chemical reactions and is a good solvent. Changing any part of the Earth system, such as the amount or type of vegetation in a region or from natural land cover to an impervious one, can affect the rest of the system. Rain and snow capture aerosols from the air. Acidic water slowly dissolves rocks, placing dissolved solids in water. Dissolved or suspended impurities determine water's chemical composition.

Current measurement programs in many areas of the world cover only a few water bodies a few times during the year. GLOBE Hydrosphere protocols will allow you to collect valuable data to help fill these gaps and improve our understanding of Earth's natural waters.



The Earth System: Energy flows and matter cycles.

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself


H. Additional resources



The Hydrosphere

GLOBE has a number of protocols that are part of the Hydrosphere Investigation. To measure the acidity of a water body, you will use one of the water pH protocols. You have the choice of using either a pH meter or pH paper in your investigations.

This slide set provides information to using a pH meter.



GLOBE Hydrosphere Measurements

Hydrosphere Study Site

Water Temperature

Water Transparency

Conductivity

pH

Mosquito Larvae

Alkalinity

Dissolved Oxygen

Salinity

Nitrates

Freshwater Macroinvertebrates

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources




What is Water pH?

pH is a measure of the relative amount of free hydrogen ions there are in the water, which determines the acidity of the water body

The concentration of the hydrogen ion $[H^+]$ activity in a solution determines the pH. Mathematically this is expressed as:

$$pH = -\log [H^+]$$



pH is reported in logarithmic units from 0-14, with 7 being neutral. Each number represents a 10x change in the acidity or alkalinity of the water.

The pH values for your water site will depend on the geology, soil and vegetation of your area as well as other inputs into your water body. Where the air masses come from that precipitate into your water body may affect the pH of the water. Most lakes and streams have pH values that range between 6.5 and 8.5. Oceans have a pH value of 8.2. Pure water not in contact with the air has a neutral pH value of 7.0.

Naturally occurring basic waters are found typically in areas where the surrounding geology is rich in minerals such as calcite or limestone.

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources




pH and Aquatic Life

pH affects most chemical and biological processes that take place in water. pH affects the solubility (amount that can be dissolved in water) and biological availability of nutrients. It also determines the degree to which potentially toxic materials, such as heavy metals, are soluble.

pH has a strong influence on what can live in the water; aquatic organisms have certain pH ranges they prefer or require. Salamanders, frogs and other amphibian life, as well as many macroinvertebrates, are particularly sensitive to extreme pH levels. Most insects, amphibians and fish are absent in water bodies with pH below 4.0 or above 10.0.

Since most organisms are sensitive to changes in water pH, scientists monitor unusual decreases or increases in the pH of water bodies. pH does not normally change a great deal, although you may find some seasonal trends due to changes in temperature, rainfall patterns, or land cover.

Importance of pH to Aquatic Life





A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.



F. Understand the data.

G. Quiz yourself

H. Additional resources

Acid pH and Aquatic Ecosystems

Without pollution, rain would have a pH of about 5.6. However, in most parts of the world rain pH is more acidic than this. Sulfur dioxide (SO_2) and nitrogen oxides (NO_x) and other compounds mix with water vapor and return to the Earth surface as acid rain. Acid rain can have a pH of



4.0 - 10.1 Limits for the most resistant fish species

Acid deposition has many harmful effects on aquatic ecosystems:

- As the pH approaches 5, invasion of non-desirable species of plankton and mosses can occur and some fish, such as smallmouth bass, begin to disappear
- Below a pH of 5, fish populations dwindle and disappear, the bottom is covered with non-decayed material, and mosses may colonize in wetlands
- Below a pH of 4.5, the water is essentially devoid of fish



A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself



H. Additional
resources

Alkaline pH and Aquatic Ecosystems

Highly alkaline (basic) waters are also hazardous for aquatic life. For example, ammonia is ten times more toxic at a pH of 8 than it is at pH 7.

- When the pH of freshwater becomes highly alkaline (e.g. 9.6), the effects on fish may include:
- death
- damage to outer surfaces like gills, eyes, and skin inability to dispose of metabolic wastes.

4.3-4.7 – nearly 10x the acidity of natural deposition.





Why Collect Water pH Data?

There are waters that are naturally more acidic when there are certain types of minerals in the water, such as sulfides. Mining can also release acid forming compounds to water bodies.

Pollution can change a water's pH, which in turn can harm animals and plants living in the water. For instance, the 2015 spill of mine waste into the Animas River, caused the Animas to have a pH of 5, that is, acidic. By using the logarithm scale, this mine-drainage water would be 100 times more acidic than neutral water.

A change in pH in a water body can be an indicator of increasing pollution or other environmental factor.



The Animas River between Silverton and Durango in Colorado, USA, within 24 hours of the 2015 Gold King Mine waste water spill Credit: Riverhugger, Wikipedia Commons⁸

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data collection!
Question 1

When determining pH, you are measuring:

- A. The relative amount of free hydrogen ions there are in the water
- B. The total dissolved solids in the water
- C. The ability of water to transmit an electrical current

What is the answer?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data collection!
Answer to Question 1

When determining pH, you are measuring:

- A. The relative amount of free hydrogen ions there are in the water 😊 **Correct!**
- B. The total dissolved solids in the water
- C. The ability of water to transmit an electrical current

Were you correct?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data collection!
Question 2

Logarithmic scale- Between pH 5 and 6, there is a ____ change in acidity.

- A. 2 x
- B. 10 x
- C. 100 x
- D. 1000 x

What is the answer?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data collection!
Answer to Question 2 !

Logarithmic scale- Between pH 5 and 6, there is a ____ change in acidity.

- A. 2 x
- B. 10 x 😊 Correct!**
- C. 100 x
- D. 1000 x

Were you correct?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data collection!
Question 3

Which of the following pH values has greater acidity?

- A. pH 4
- B. pH 7
- C. pH 9

What is the answer?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data collection!
Answer to Question 3

Which of the following pH values has greater acidity?

A. pH 4 - 😊 Correct!

B. pH 7

C. pH 9

Were you correct?

Let's now look at the data collection procedures!

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



pH Protocol using pH Paper: What do you need to start?

When	Suggested Frequency: Weekly
Where	Hydrosphere Study Site
Time needed	10 minutes
Prerequisites	Hydrosphere Study Site Definition
Key instruments	pH paper, Electrical Conductivity meter
Skill level	all
Necessary Documents	<ul style="list-style-type: none">• <u>Hydrosphere Investigation Data Sheet</u>• <u>pH Water Protocol using pH Paper:</u>• <u>Electrical Conductivity Greater than 200 mS/</u>• <u>Using pH Paper (electrical conductivity less than 200 mS/cm)</u>

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Overview of Water pH Protocol

A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources

- pH of a water body can be measured using either a pH meter or pH paper. The accuracy of either method depends on the **electrical conductivity** of the water. The electrical conductivity needs to be at least 200 $\mu\text{S}/\text{cm}$ for these methods to report accurately.
- If you are sampling Ocean or brackish water, you can assume that the electrical conductivity of your sample is greater than 200 $\mu\text{S}/\text{cm}$. If you are not sure if the fresh water at your Hydrosphere Study Site has a conductivity value high enough for the measurement technique (paper or meter), you will need to measure the **electrical conductivity** before taking your pH measurements. After you know the electrical conductivity value of the water, use the appropriate pH field guide.
- For more information, see the [Electrical Conductivity Field Guide Protocol](#)



Simultaneous or Prior Investigations Required to do Water pH Measurements: Hydrosphere Study Site Description

You will need to define your **Hydrosphere Study Site**. A **Hydrosphere Study Site** can be any surface water site that can be safely visited, although natural waters are preferred.

The **Hydrosphere Investigation Data Sheet** is used to record all the hydrosphere measurements, including Water Transparency. You will also want to map your Hydrosphere Site at some point.

To define your study site you will need these documents:

- [Selecting and Documenting your Hydrosphere Study Site](#)
- [Hydrosphere Investigation Data Sheet](#)
- [Mapping your Hydrosphere Study Site Field Guide](#)

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Site Selection and Sampling- Hydrosphere Study Site

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources

If the selected study site is a moving body of water (i.e. stream or river), locate your sampling site at a riffle area as opposed to still water or rapids. This will provide a more representative measurement of the water in the stream or river. If the selected study site is a still body of water i.e. a lake or reservoir), find a sampling site near the outlet area or along the middle of the water body.



Avoid inlet areas. A bridge or a pier are good choices. If your water body is brackish or salty, you will need to know the times of high and low tide at a location as close as possible to your study site.



Prior Investigations Required: Electrical Conductivity

To obtain accurate pH measurement of your water sample, the electrical conductivity of your sample must be at least 200 $\mu\text{S}/\text{cm}$. After completing the Electrical Conductivity Protocol, refer to your results and select the pH protocol that is appropriate. In this tutorial, we will review all of the protocols, so you can see the differences in procedure between them. You will see that if the Electrical Conductivity is less than 200 $\mu\text{S}/\text{cm}$, you will need to do an extra step, and add salt to the sample until it is at least 200 $\mu\text{S}/\text{cm}$. This is true for both pH paper and pH meter methods.

This slide set includes these protocols:

- [I. Using pH Paper \(electrical conductivity greater than 200 \$\mu\text{S}/\text{cm}\$ \)](#)
- [II. Using pH Paper \(electrical conductivity less than 200 \$\text{mS}/\text{cm}\$ \)](#)

A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources



Start your Field and Lab work with Safety Steps

A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources

Safety is important when conducting the Hydrosphere protocols. While you will need to use your judgment in selecting only hydrosphere study sites that are safe to access and sample, additional precautions are needed.

Students should wear protective gloves and goggles when handling water samples and chemicals to avoid danger from splashes.

When doing GLOBE Hydrosphere Protocols, it is important to protect students from exposure to biting mosquitoes. Ask your students to wear clothes that cover the body so there is little bite area exposed. It is also advisable to apply insect repellent if you are sampling during the mosquito breeding season.



SAFETY be sure students wear gloves and goggles during your investigations





A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources

I. pH Water Protocol using pH Paper: Electrical Conductivity Greater than 200 $\mu\text{S}/\text{cm}$



I. pH Water Protocol using pH Paper: Electrical Conductivity Greater than 200 $\mu\text{S}/\text{cm}$ (Slide 1/3)

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources

1. Fill in the top part of the Hydrosphere Investigation Sheet
2. In the pH section of the Data Sheet, check the box next to "pH paper"
3. Put on protective gloves and goggles

Hydrosphere Investigation

Data Sheet

School name: _____ Class or group name: _____

Name(s) of Student(s) collecting data: _____

Measurement Time: *

Year: _____ Month: _____ Day: _____ Time: _____ (UT) Time: _____ (Local)

Name of Site: _____

Water State: (check one) *

☐ Normal ☐ Flooded ☐ Dry ☐ Frozen ☐ Unreachable

Note: If Normal is selected, continue below; all other selections stop here

Sky Conditions (Check one):

- ☐ Clear (no Clouds Visible)
☐ Clouds Visible (1% to 100% Covered by Clouds or Contrails)
☐ Obscured (More than 25% of the Sky is not Visible)

Note: selecting **Obscured** will prevent data entry on clouds and contrails; therefore skip the cloud type and cover and the contrail type and cover sections and proceed to the Obscured section. If clouds and contrails are visible in non-obscured areas of the sky, these data can be entered in the Metadata field.

Water pH: Measured with: (check one) ☐ pH Paper ☐ pH Meter

If salt added, conductivity ($\mu\text{S}/\text{cm}$)	pH
1.	
2.	
3.	

Value of buffers used: ☐ pH 4 ☐ pH 7 ☐ pH 10 (Check all used)

Comments: _____



I. pH Water Protocol using pH Paper: Electrical Conductivity Greater than 200 $\mu\text{S}/\text{cm}$ (2/3)

4. Rinse the beaker with sample water **three** times

5. Fill the beaker halfway with your water sample

6. Follow instructions that come with your pH paper for testing the sample

7. Record your pH on the Data Sheet as Observer 1



A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



I. pH Water Protocol using pH Paper: Electrical Conductivity Greater than 200 $\mu\text{S}/\text{cm}$ (3/3)

5. **Repeat** steps 4-6 two more times, using fresh sample water and new pH paper.
6. Record your data, and find the average of the three observations:
 - **Calculate Average:**
 - **Observer 1 + Observer 2 + Observer 3**
10. Discard used pH paper and gloves in waste container. Rinse beaker with distilled water using the wash bottle.
11. Enter your data on the GLOBE Website.



***End of data collection for this Water pH Protocol ***



Check to make sure that each observation is within 1.0 pH units of the average. If they are not within 1.0 units of the average, repeat the measurements

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



- A. What is water pH?
- B. Why collect water pH data?
- C. How your measurements can help
- D. How to collect your data.
- E. Entering data on GLOBE Website.
- F. Understand the data.
- G. Quiz yourself
- H. Additional resources

II. pH Water Protocol using pH Paper: Electrical Conductivity Less than 200 $\mu\text{S}/\text{cm}$



II. pH Water Protocol using pH Paper: Electrical Conductivity Less than 200 $\mu\text{S}/\text{cm}$



- Assemble Necessary Additional Documents:
- [pH Water Protocol using pH Paper: Electrical Conductivity Less than 200 \$\mu\text{S}/\text{cm}\$](#)



*Check to make sure that each observation is within 1.0 pH units of the average.
If they are not within 1.0 units of the average, repeat the measurements*

A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources



A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources

II. pH Water Protocol using pH Paper: Electrical Conductivity Less than 200 $\mu\text{S}/\text{cm}$ (slide 1/4)

1. Fill in the top part of the Hydrosphere Investigation Sheet. In the pH section of the Data Sheet, check the box next to “pH paper”.

2. Put on protective gloves and goggles

3. Rinse tweezers in sample water and dry with paper towel

4. Rinse the beakers with sample water three times.





A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources

II. pH Water Protocol using pH Paper: Electrical Conductivity Less than 200 $\mu\text{S}/\text{cm}$ (slide 1/4)

5. Fill one beaker or cup with about 50 mL of sample water

6. Using the tweezers, place one crystal of salt in the sample water. (If you do not have salt crystals, use several grains of table salt (a pile about 2mm wide at the bottom))

7. Stir thoroughly with stirring rod or spoon.



*A note regarding salt crystals. Crystal of about 0.5 – 2.0 mm in diameter are much easier to work with than the very finely ground “table salt” used in some countries. In North America, the larger salt crystals are often marketed as “sea salt”



II. pH Water Protocol using pH Paper: Electrical Conductivity Less than 200 $\mu\text{S}/\text{cm}$ (slide 2/4)

- 8. Measure the electrical conductivity of the treated sample water (with the added salt) using the **Electrical Conductivity Protocol**
 - a. If the electrical conductivity is at least 200 $\mu\text{S}/\text{cm}$, record value on
 - Data Sheet. Go to **step 9**.
 - b. If the electrical conductivity is still less than 200 $\mu\text{S}/\text{cm}$, go back to **step 6** and repeat until you get a value that is at least 200 mS/cm. Record conductivity value on Data Sheet.



A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



II. pH Water Protocol using pH Paper: Electrical Conductivity Less than 200 $\mu\text{S}/\text{cm}$ (slide 3/4)

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources

9. Follow instructions that come with your paper for testing of the pH of the sample.

10. Record your pH on the Data Sheet as Observer 1.

11. Repeat steps 3-9 using new water samples and new pieces of paper. Record the the data on the Data Sheet as Observer 2 and Observer 3.

12. Find the average of the three observations.

Calculate Average: Observer 1 + Observer 2 + Observer 3

3



If you are unable to get readings within 1 pH unit in your three samples, there may be problems with your pH paper.





II. pH Water Protocol using pH Paper: Electrical Conductivity Less than 200 $\mu\text{S}/\text{cm}$ (slide 4/4)

13. Check to make sure that each observation is within 1.0 pH units of the average. If they are not within 1.0 units of the average, repeat the measurements.

14. Discard used pH paper and gloves in a waste container. Rinse the beaker with distilled water.

15. Enter your data on the GLOBE Website.

End of data collection for this pH Water Protocol



A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data reporting!
Question 4

Which of the following conditions must be met to ensure the pH measurement is accurate?

- A. An electrical conductivity of 200 $\mu\text{S}/\text{cm}$ or greater
- B. A temperature of greater than 30 degrees C
- C. Both A and B
- D. None of the above, you can just stick the pH paper in the water and take the measurement, and it will be accurate

What is the answer?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data reporting!
Answer to Question 4

Which of the following conditions must be met to ensure the pH measurement is accurate?

- A. An electrical conductivity of 200 $\mu\text{S}/\text{cm}$ or greater-**correct** 😊**
- B. A temperature of greater than 30 degrees C
- C. Both A and B
- D. None of the above, you can just stick the pH paper in the water and take the measurement, and it will be accurate

Were you correct?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data reporting!
Question 5

How many replicate measurements should be taken for each sample you are testing for pH to ensure accuracy?

- A. 2
- B. 3
- C. 4-6

What is the answer?

- A. What is water pH?
- B. Why collect water pH data?
- C. How your measurements can help
- D. How to collect your data.
- E. Entering data on GLOBE Website.
- F. Understand the data.
- G. Quiz yourself
- H. Additional resources



Let's do a quick review before moving onto data reporting!
Answer to Question 5

How many replicate measurements should be taken for each sample you are testing for pH to ensure accuracy?

- A. 2
- B. 3- 😊 Correct!**
- C. 4-6

Were you correct?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data reporting!
Question 6

To accept your measurements for reporting to GLOBE, Each of the independent pH measurements should be

- A. within 1.0 pH units of the average of the replicate measurements
- B. within $-\log_{10}$ of the average
- C. within 1.0 pH units of pure water

What is the answer?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Let's do a quick review before moving onto data reporting!
Answer to Question 6

To accept your measurements for reporting to GLOBE, Each of the independent pH measurements should be

- A. within 1.0 pH units of the average of the replicate measurements - 😊 Correct!**
- B. within $-\log_{10}$ of the average
- C. within 1.0 pH units of pure water

Were you correct?

Let's move on to GLOBE data entry and visualization!

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Submit your Data to GLOBE

A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

**E. Entering
data on GLOBE
Website.**

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources

- [Live Data Entry](#): Upload your data to the official
- GLOBE science database
- Email Data Entry: Send data in the body of your email (not as an attachment) to DATA@GLOBE.GOV
- Mobile Data App: Download the GLOBE Science Data Entry app to your mobile device and select the right option.
- **For Android** via [Google Play](#)
- **For IOS** via the [App Store](#)





Entering your data via Live Data Entry or Data Entry Mobile App- Step 1

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources

Identify your
Sampling site

Select
"Integrated Hydrology"
and "New observation"

Welcome to the GLOBE data entry site. [X]

My Bookmarks

You have not bookmarked any investigations yet. Expand the organizations and click the stars next to the investigations to create a bookmark.

My Organizations and Sites

- [University of Nebraska-Lincoln GLOBE v-School](#) [Add site]

+ [Lefthand Creek](#)
Latitude 40, Longitude -105, Elevation 1600m
[Edit site] [Delete site]

Hydrology

Freshwater Macroinvertebrates ★
[New observation] [Past observations]

Integrated Hydrology ★
[New observation] [Past observations]

Mosquitoes ★
[New observation] [Past observations]

Annotations: An arrow points from "Identify your Sampling site" to the "Lefthand Creek" site entry. Two arrows point from "Select 'Integrated Hydrology' and 'New observation'" to the "Integrated Hydrology" section and its "New observation" button respectively.



Entering your data via Live Data Entry or Data Entry Mobile App- Step 2

1. Select
water body
state

2. Select
protocol

**Be sure to enter data
as pH paper**

3. Enter each
measurement
and click "add"

4. Click to
send data

**You are done! Want to
check who else has
submitted pH data
using the GLOBE
Visualization System?**



A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources

Visualize and Retrieve Water pH Data- Step 1

GLOBE provides the ability to view and interact with data measured across the world. Select our [visualization tool](#) to map, graph, filter and export pH data that have been measured across GLOBE protocols since 1995. Here are screenshots steps you will use when you use the visualization tool:



[Link](#) to step-by-step tutorial on using the GLOBE Data Visualization System



Visualize and Retrieve Water pH Data- Step 2

Select the date for which you need pH data, add layer and you can see where data is available.



Locations where pH
data is available for
the week you
selected

A. What is
water pH?

B. Why collect
water pH data?

C. How your
measurements
can help

D. How to
collect your
data.

E. Entering
data on GLOBE
Website.

F. Understand
the data.

G. Quiz
yourself

H. Additional
resources



Visualize and Retrieve Water pH Data- Step 3

Select the sampling site for which you need pH data, and a box will open with data summary for that site.



Clicking on a location will open to a map note providing pH data for that location and time. Follow instructions in the tutorial to download data as a .csv file for analysis.



Review questions to help you prepare to conduct the Hydrosphere pH Protocol

1. What is the importance of pH to aquatic life?
2. What is a logarithmic scale? Why is it a useful way to report pH?
3. True/False pH affects the solubility and biological availability of nutrients.
4. In a water body, what happens to aquatic life in waters with pH values below 4.0 or above 10.0?
5. What other measurement, in addition to pH, must you do to ensure that your pH paper or meter is reporting accurately?
6. What are the safety precautions you should take when doing any of the hydrology protocols?
7. What is the acceptable range of error of the three replicate samples you take?
8. What kinds of environmental events could change water pH?
9. Which pH value is more acid: a pH value of 1 or a pH value of 14?
10. What is the pH value of pure water?

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



Are you ready to take the quiz?

- A. What is water pH?
- B. Why collect water pH data?
- C. How your measurements can help
- D. How to collect your data.
- E. Entering data on GLOBE Website.
- F. Understand the data.
- G. Quiz yourself
- H. Additional resources

- You have now completed the slide stack. If you are ready to take the quiz, sign on and take the quiz corresponding to **Water pH Protocol**.
- You can also review the slide stack, post questions on the Hydrosphere discussion area, or look at the FAQs on the next page.
- When you pass the quiz, you are ready to take **Water pH Protocol** measurements!



FAQ: Frequently Asked Questions

Does water temperature affect my pH reading?

A change in water temperature can actually change the pH value of your water. Since we want to know the actual pH value, we do not correct for this change.

Temperature can also affect the performance of the meter. The electrode is designed so there is no temperature sensitivity when the pH is 7.0. As the pH moves away from this value, the water temperature affects meter accuracy. Meters with automatic temperature compensation (ATC) correct for the temperature of the water at values above and below 7.0 by a factor of 0.003 pH/°C/pH unit away from pH 7. They correct for meter error.

Does high salt concentration affect pH?

Salt concentration can affect pH. As salt concentration increases, pH can increase. This is not a linear relationship, but can be important in estuaries, where the salinity varies with the tide. Taking into account salinity or conductivity data may be useful in understanding variations in your pH measurements.

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



FAQ: Frequently Asked Questions-2

- **Why may pH measurements be inaccurate in low conductivity waters?**
 - To measure the hydrogen ion concentration, you are actually measuring the potential of the hydrogen ions. Other ions have to be present to pass the current to make this measurement. When they are at too low of a concentration the meter slowly drifts and if the drift is really slow, the meter locks in on an incorrect measurement
- **Can I use a pH meter that connects to my Smart Phone?**
 - Yes, pH meters that connect to iPhones, iPads and other Smart devices can be considered pH meters. For most of these meters an app is required.

A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources



A. What is water pH?

B. Why collect water pH data?

C. How your measurements can help

D. How to collect your data.

E. Entering data on GLOBE Website.

F. Understand the data.

G. Quiz yourself

H. Additional resources

Please provide us with feedback about this module. This is a community project and we welcome your comments, suggestions and edits! Comment here: [eTraining Feedback](#)
Questions about module content? Contact GLOBE eTraining: rlow@ucar.edu

Credits

Slides:

Russanne Low, Ph.D., University of Nebraska-Lincoln, USA

Rebecca Boger, Ph.D., Brooklyn College, NYC, USA

Photos: Russanne Low

Art: Jenn Glaser, *ScribeArts*

More Information:

[The GLOBE Program](#)

[NASA Wavelength](#) NASA's Digital Library for Earth and Space Education

[NASA Global Climate Change: Vital Signs of the Planet](#)

The GLOBE Program is sponsored by these organizations:



Version 12/1/16. If you edit and modify this slide set for use for educational purposes, please note "modified by (and your name and date)" on this page. Thank you.